

Microbiological Drinking Water Analysis					
Test	Method	Bottle	Perservative	Hold Time	Turn Time
Total Coliform MMO-MUG	SM 9223-18	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
E-Coli MMO-MUG	SM9223 B	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
Total Coliform Quanti Tray	SM9223 B-QT	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
Total Coliform Membrane filter	SM 9222 B G	100ml Sterile	Sodium Thiosulfate	30 hours	96 hours
Fecal Coliform Membrane filter	SM 9222 D	100ml Sterile	Sodium Thiosulfate	30 hours	36 hours
E.Coli Membrane filter	EPA 1603	100ml Sterile	Sodium Thiosulfate	30 hours	96 hours
Total Plate Count	SM 9215 B	100ml Sterile	Sodium Thiosulfate	24 hours	96 hours